

**PROFESSIONAL  
PROCEDURE  
GUIDE**



CLIA-WAIVED A1C TESTING

METRIKA

METRIKA


CLIA-WAIVED A1C TESTING



**PROFESSIONAL  
PROCEDURE  
GUIDE**

**BEFORE YOU BEGIN**

**Make sure all parts are the same lot number.** Always run the test with all parts of the test kit at room temperature (18°–28°C, 64°– 82°F). If the kit has been recently at high temperatures (above 82°F) or in the refrigerator allow **all** parts to come to room temperature in their sealed foil pouches for at least one hour before running the test.



Avoid running the test in direct sunlight, on hot or cold surfaces, or near sources of heat or cold.


**PREPARE**

Prepare the following items. Have them all at room temperature:

**1. The Monitor**



**2. Sample Dilution Kit (Open Pouch)**





**3. Test Cartridge (Wait to open)**



**GET BLOOD**

**OPEN POUCH #1**, then squeeze the folded Tube Holder to open it. Remove cap. Place Tube in Holder.

**CLEAN AND DRY** the patient's finger. Lance the finger to obtain a large (10 µL) drop of blood. Do not milk the finger.




Do not underfill, overfill, or leave a hanging drop.

**NOTE:** For Venous Draw Procedure: Obtain a blood sample by standard venipuncture technique in an EDTA (purple top) tube. Mix the sample well before testing. A standard laboratory precision pipet may be used to transfer 10 µL from an EDTA tube into the Dilution Buffer tube, instead of the capillary pipet.


**ADD BLOOD TO TUBE**

**PUT CAPILLARY TIP** into the liquid in the Tube. Leaving it submerged, squeeze the bulb firmly two or three times to rinse all of the blood from the Capillary into the Tube.



**RECAP TUBE, SHAKE!**

**RECAP THE TUBE** tightly and shake it vigorously 6-8 times. It's OK to have some bubbles. The diluted sample will be red-orange in color. Replace the Tube in the Tube Holder.



**INSERT CARTRIDGE**

**OPEN POUCH #2** and insert Cartridge into Monitor immediately. Place on a flat surface. The display will turn on and show "WAIT" until the monitor performs internal QC checks. After these QC checks are complete, the monitor will show "SMPL".

The Monitor is now ready for use. (If you see any other message on the display, go to "Troubleshooting").



Match codes

**NOTE:** The code on the Monitor must match the code on the Cartridge. If the code numbers do not match, DO NOT continue with the test.

**DILUTED SAMPLE MUST BE ADDED WITHIN 2 MINUTES.**

**FILL DROPPER**

1. Squeeze bulb
2. Lower into vial
3. Release
4. Rapidly remove from vial. DO NOT SQUEEZE
5. Check for good fill

**ADD DILUTED SAMPLE**

**REST** dropper in port



Squeeze bulb fully and quickly



Pull away, while still squeezing bulb



**In summary: Add a full barrel of diluted sample quickly to the sample port.**

**READ RESULT**

The display will start counting down from 5 minutes, then show a test result alternating with the letters "QC OK". For any other error messages, please see list of error codes on reverse side. Do not handle the Monitor until the result is displayed. The result remains displayed for 60 minutes or until the next Test Cartridge is inserted.



**REMOVE AND DISCARD CARTRIDGE  
KEEP MONITOR**

**SAVE MONITOR**

**THE MONITOR IS REUSABLE** until all the Cartridges are used up. To run a new test, start with a new Sample Dilution Kit and Test Cartridge and go to "Prepare".



The monitor will display 3 TL, 2 TL, 1 TL and 0 TL, alternating with "QC OK" and the result when there are 3, 2, 1 and 0 tests left, respectively. After all Cartridges in the kit are used, the Monitor will expire. (If you insert a new Cartridge, the display will show "0 TL" (Zero Tests Left) for five minutes and then shut down.)

**Metrika, Inc.**  
510 Oakmead Parkway  
Sunnyvale, CA 94085-4022  
tel 877 212 4968  
fax 408 524 6595  
www.metrika.com

# A1CNow+ PROFESSIONAL-USE PRODUCT INSERT

## Intended Use

The A1cNow+® test provides quantitative measurement of the percent of glycated hemoglobin (%A1C) levels in capillary (fingerstick) or venous whole blood samples. The test is for professional use to monitor glycemic control in people with diabetes.

## Summary and Explanation

High levels of blood glucose result in over-glycation of proteins throughout the body including hemoglobin.<sup>1</sup> Glycation of hemoglobin can occur at the amino termini of the alpha and beta chains, as well as other sites with free amino groups<sup>1</sup>. Hemoglobin A under goes a slow glycation with glucose that is dependent on the time-average concentration of glucose over the 120-day life span of red blood cells.

The most prevalent and well-characterized species of glycated hemoglobin A is A1C, making up approximately 3% to 6% of total hemoglobin in healthy individuals.<sup>1</sup> The correlation of A1C and blood glucose levels make it a useful method of monitoring long-term blood glucose levels in people with diabetes.<sup>2</sup> Previous studies, such as the Diabetes Complications and Control Trial (DCCT) and the United Kingdom Prospective Diabetes Study (UKPDS), used glycated hemoglobin as a way to measure overall glycemic control during the studies. These studies, and others, have shown that tight glycemic control is associated with fewer diabetes-related complications (e.g., vision problems, cardiovascular problems, and kidney problems).<sup>3</sup> The National Glycohemoglobin Standardization Program (NGSP) was established to assure traceability of hemoglobin A1C (A1C) results to the DCCT. Studies show a direct relationship from %A1C to average blood glucose (MBG) levels. For every 1% change in A1C there is a change of about 30 mg/dl in MBG.<sup>4</sup> The formula used to calculate the mean (average) blood glucose levels from the A1C levels is MBG = (31.7 x HbA1c) - 66.1. To convert to mean plasma glucose (MPG) use<sup>5</sup> MPG = MBG x 1.11.

A1C can be measured by a variety of techniques, and over the past decade they have expanded to include point-of-care assays. Point-of-care assays are well suited to environments such as physicians' offices and clinics, because they are generally easy to perform, require no laboratory equipment, and provide rapid turn-around-time from sampling to result.<sup>6</sup> This immediate feedback of results enhances physician/patient interaction and, therefore better enables disease management.<sup>7</sup>

## Principle of the Assay

Metrika has developed an enabling technology called MODM™ (Micro-Optical Detection Method) that incorporates microelectronics, optics, and dry-reagent chemistry strips within a reusable, self-contained, integrated hand held monitor and a single-use test cartridge. An unmeasured whole blood mixture (diluted) is directly applied to the sample port, and results are displayed in numeric form on the monitor's liquid crystal display after 5 minutes. Having no switches or buttons, the monitor self-activates upon insertion of the Test Cartridge.

The A1cNow+ monitor utilizes both immunoassay and chemistry technology to measure A1C and total hemoglobin, respectively. Upon the addition of a diluted blood sample, blue microparticles conjugated to anti-A1C antibodies migrate along the reagent strips. The amount of blue microparticles captured on the strips reflects the amount of A1C in the sample.

For the total hemoglobin (Hb) portion of the test, the sample diluent converts Hb to met-Hb. The intensity of met-Hb color measured on the reagent strips is proportional to the concentration of hemoglobin in the sample. Test results are expressed as %A1C (A1C ÷ total Hb x 100).

Calibration of the A1cNow+ is performed with a set of blood samples that have been value-assigned by a National Glycohemoglobin Standardization Program (NGSP) certified laboratory using an NGSP reference method. Total Hb calibration values for those samples are obtained with a Total Hb analyzer (HemoCue Hemoglobin Test System, HemoCue, Inc., Lake Forest, CA). The calibration of the A1cNow+ test is thus traceable to the NGSP and to an NGSP Certified Network reference method.

## Specimen Collection and Storage

**Note: No fasting or special diet is necessary**

## Fingerstick

The A1cNow+ test requires 10 microliters (µL) of whole blood (1 large drop). Fingerstick blood is obtained by standard techniques with any lancing system. If alcohol is used for cleansing, be sure the finger is completely dry before lancing.

## Venipuncture

Venous blood is to be collected in an EDTA tube ("Purple Top"). Blood should be well-mixed, and tested at room temperature. Venous blood samples are stable for up to 8 hours at room temperature, and up to 14 days if refrigerated (2-8°C).

## Warnings and Precautions

- For in vitro diagnostic use only.
- Carefully read and follow the Professional Procedure Guide (on reverse) to ensure proper test performance.
- If refrigerated, bring sealed pouches and monitor to room temperature for one hour.
- The A1cNow+ Monitor and Test Cartridges should not be used if either are cracked or broken.
- The Test Cartridges should not be used if the foil pouch is damaged.
- Add sample to A1cNow+ Test Cartridge within 2 minutes after pouch is opened.
- All components of the A1cNow+ system are potentially bio-hazardous. Dispose of as biohazardous waste.
- The Dilution Buffer contains ferricyanide in a buffered detergent solution. Do Not Ingest. In case of contact with skin or eyes, flush the area with large amounts of water.
- Do not reuse Test Cartridges or Sample Dilution Kits.

**Do not mix Monitors with Cartridges & Sample Dilution Kits from different lots.**

## Kit Storage and Stability

• Pouched Test Cartridges, A1cNow+ Monitors, and Sample Dilution Kits may be stored at room temperature (18-28°C) for up to **three months** prior to use. Monitors, Test Cartridges, and Dilution Kits must be thrown away if not used within the **three months**.

• The Monitors, Test Cartridges, and Sample Dilution Kits may be used until the expiration date printed on the box and pouches when stored refrigerated (2-8°C). Monitors, Test Cartridges, and Sample Dilution Kits must be thrown away if not used by the expiration date.

• Leave all components in their sealed pouches until use. If refrigerated, ensure pouches are at room temperature before use.

• Do not mix pouches and Monitors from different lots.

## Package Components

- A1cNow+ Monitor (1)
- A1cNow+ Test Cartridges (10, or 20) Each Test Cartridge includes the following chemistries: antibody to HbA1c, antigen conjugate that binds to the antibody, and membranes.
- Sample Dilution Kit (10, or 20), each containing:
  - Tube (1), containing 0.69 mL of buffered detergent solution with ferricyanide
  - Capillary (1)
  - Dropper (1)
  - Tube Holder (1)
- Product insert (1)
- Patient result labels (10, or 20)

## Materials Required but Not Supplied

- Fingerstick sample: lancet, or other blood fingerstick collection device or,
- Venous sample: EDTA tube ("purple top"), venous collection supplies,
- Gauze pad or cotton ball
- Bandage

## Result Interpretation

Percent A1C monitors glucose control over the last three months. About 50% of the A1C result is from the past 30 days; about 25% is from the past 30-60 days and about 25% is from the past 60-120 days.<sup>4</sup> Depending on the test methodology used, laboratory methods show that the reference range of the A1C test is approximately 4.0-6.5% A1C, and 6% to 9% in people with well to moderately controlled diabetes.<sup>1</sup> Levels can be as high as 20% in people with poorly controlled diabetes.<sup>8</sup> The American Diabetes Association's (ADA's) most recent Clinical Practice Recommendation for diabetes specifies a treatment goal for patients in general of less than 7% with a treatment goal for the individual patient of as close to normal (less than 6%) as possible without significant hypoglycemia.<sup>9</sup>

## Troubleshooting

See the table below for a description of A1cNow+ operating and error codes (OR =Out of Range; QC = Quality Control, E= Monitor Error)

MESSAGE	DESCRIPTION AND RESOLUTION
<b>OR 1</b>	The blood sample may have too little hemoglobin (less than 20% hematocrit), or there was under-sampling of whole blood.* You may wish to check hemocrit by another method.
<b>OR 2</b>	The blood sample may have too much hemoglobin (greater than 60% hemocrit), or there was over-sampling of whole blood.* You may wish to check hemocrit by another method.
<b>OR 3</b>	The blood sample may have too little A1C, or there was under-sampling of whole blood.*
<b>OR 4</b>	The blood sample may have too much A1C, or there was over-sampling of whole blood.*
<b>OR 5</b>	The monitor temperature is below 18°C (64°F). Repeat the test at room temperature (18-28°C).
<b>OR 6</b>	The monitor temperature is above 28°C (82°F). Repeat the test at room temperature (18-28°C).
<b>&lt;4.0</b>	The %A1C is less than 4%.
<b>&gt;13.0</b>	The %A1C is greater than 13%.
<b>QC 2</b>	Occurs when you insert a cartridge that already has sample added to it. Do not remove and reinsert a cartridge after adding sample.

<b>QC 6</b>	Sample was added to cartridge before "SMPL" display. This counts down one test on the monitor. Remove and discard cartridge. To avoid this error, do not add sample until the "WAIT" prompt clears and "SMPL" appears.
<b>QC 7</b>	The cartridge remained in the monitor without sample addition for 2 minutes after "SMPL" prompt. This counts down one test on the monitor. Discard the cartridge and insert a fresh one when you are ready to add the diluted sample.

<b>QC 50 to 51</b> <b>QC 55 to 56</b>	Insufficient sample was added to the cartridge or it was added too slowly. Remove and discard cartridge. To avoid this error, add a full dropper of diluted sample to the sample port in less than a second.
<b>All other QC codes</b>	The quality control checks did not pass. Call Metrika Technical Support toll free at 877-METRIKA x522. The test will have to be repeated with another Test Cartridge and Sample Dilution Kit.

<b>E1 to E99</b>	The Monitor has a Fatal Error. Call Metrika Technical Support toll-free at 877-METRIKA x522.
------------------	--

*Carefully repeat the test using a new Test Cartridge and a new Sample Dilution Kit.
--

This test is WAIVED under the Clinical Laboratory Improvement Amendments of 1988 (CLIA).

If a laboratory modifies the test instructions, the test will no longer be considered waived.

## Limitations

- This test is NOT for the screening or diagnosis of diabetes.
- If the patient has high levels of Hemoglobin F, Hemoglobin S, Hemoglobin C, or other hemoglobin variants, the A1cNow+ system may report incorrect results.
- Any cause of shortened red cell survival (e.g., hemolytic anemia or other hemolytic diseases, pregnancy, recent significant blood loss, etc.) will reduce exposure of red cells to glucose. This results in a decrease in %A1C values. Percent A1C results are not reliable in patients with chronic blood loss and consequent variable erythrocyte life span.
- This test is designed to be run at 18-28°C (64-82°F) and 15-80% humidity. Using the monitor outside this temperature range will give an error code.
- This test is not a substitute for regular doctor visits and blood glucose monitoring.
- As with any laboratory procedure, a large discrepancy between clinical impression and test results usually warrants investigation.

## Controls

Each A1cNow+ monitor performs over 25 internal chemical and electronic quality control checks, including potential hardware and software errors (e.g., cartridge alignment, programming), and potential reagent strip errors (e.g., insufficient sample volume, invalid calculations). The monitor has been programmed to report an error code if these quality checks are not passed.

If external quality control testing is desired, commercial controls may be purchased from other vendors. Please contact Metrika Customer Service for recommendations. Metrika recommends that external controls be tested at the following times:

- Whenever laboratory or room conditions have been above 28°C if stored at room temperature.
- To perform training or retraining of testing personnel.
- Whenever A1cNow+ results do not match other clinical findings or symptoms.

Good laboratory practices include a complete quality control program. This entails proper sample collection and handling practices, ongoing training of testing personnel, ongoing evaluation of control results, proper storage of test kits, etc. A permanent record of control results should be retained.

## Performance

### Expected Values (non-diabetic population)

The expected normal range for %A1C using the A1cNow+ system was determined by testing blood samples from 118 presumptively non-diabetic individuals (fasting glucose levels <127 mg/dL) across three US sites. The population included 33 males and 85 females, and an age range from 19 to 76 with a mean age of 43. The mean %A1C result was 5.2% ±0.71% (1 SD). The 95% confidence limits were 3.9% to 6.5%. These values are similar to those reported in the literature. Each laboratory should determine its own reference range to conform to the population being tested.

## Linearity

Studies were performed to evaluate the linearity of the A1cNow+ system across its dynamic range. Clinical samples representing low and high %A1C levels were identified, and were mixed in various proportions into nine preparations. These samples were tested in replicates of at least five (n = 5). The observed results were compared to the expected results and analyzed in terms of percent recovery. The test is linear for %A1C levels between 4% and 13%, and produces reliable results with hematocrits between 20% and 60% packed cell volume (PCV).

### Interference Testing/Specificity

Studies were performed to assess the effect of common test interferents, various common over-the-counter therapeutic agents, and oral antihyperglycemic agents commonly used to treat Type II diabetes. Two levels of %A1C (low and high, approximately 4% and 10%, respectively) were tested. See table below.

INTERFERENT	TEST CONCENTRATION
Bilirubin (unconjugated)	20 mg/dL
Triglyceride	3000 mg/dL
Hemoglobin	500 mg/dL
Acetaminophen	80 mg/dL
Ascorbic acid	5mg/dL
Ibuprofen	120 mg/dL
Acetylsalicylic acid	1mg/dL
Glyburide (glibenclamide)	240 mg/dL
Metformin (1,1-dimethylbiguanide HCl)	25 mg/dL

The studies showed no effect from any of these potential interferents at concentrations up to approximately 5-times their normal levels or therapeutic doses.

Studies showed no interference from modified hemoglobins, including labile glycated hemoglobin when tested at two levels of %A1C (low and high, approximately 5% and 11% respectively). The modified hemoglobins, and the levels evaluated, were: labile hemoglobin with 1400 mg/dL glucose, carbamylated hemoglobin at a final concentration of 5 mM potassium cyanate, and acetylated hemoglobin at a final concentration of 14 mM acetylsalicylic acid.

There were mixed results from the testing of high levels of Hemoglobin F, Hemoglobin S, and Hemoglobin C. Unreliable results may be obtained from patients with elevated levels of variant hemoglobins.

## Precision

Precision testing was done under a specialized protocol. Following this protocol, two whole blood samples, one of approximately 6 %A1C (low), and one of approximately 9 %A1C (high), were tested over 20 days and four runs per day, for a total of 80 assays per level. The overall imprecision (including within-day and between-day) was 3.00% CV at the low level and 4.02% CV at the high level. This performance meets the requirements of NGSP certification.

## Accuracy

Accuracy studies were conducted with 189 diabetic and non-diabetic subjects across three US sites. Fingerstick sampling was performed on each subject for testing with A1cNow+, and venous blood was collected from each subject for comparative testing using an NGSP-certified method. A1cNow+ results were compared to the NGSP reference

results. The A1C results ranged from 5.0 %A1C to 12.8 %A1C, with a mean of 7.3 %A1C (reference results). Data analysis consisted of least squares linear regression (x = reference results), bias calculation, and Bland Altman limits. The data are provided below.

## A1cNow+ Fingerstick Comparative Testing

(NGSP-certified method is the Tosoh A1c 2.2 Plus)

n	189	Bias at 6% A1C (% difference)	5.89 (- 1.83%)
Slope	1.02	Bias at 7% A1C (% difference)	6.91 (+0.14%)
y-intercept	- 0.23	Bias at 9% A1C (% difference)	8.95 (- 0.56%)
"r"	0.95	Avg. % diff.	- 1.23%

The results showed that the accuracy of A1cNow+, with fingerstick samples, was, on average, 99%. This means that, on average, a true 7 %A1C could read approximately 6.9 %A1C. An individual A1cNow+ result may differ by as much as -1.0 %A1C to +0.8 %A1C from the true result. This represents the 95% confidence limits of a Bland-Altman plot.

## A1cNow+ System Venous Comparative Testing

(NGSP-certified method is the Tosoh A1c 2.2 Plus)

In a separate study, venous blood was collected from 50 diabetic subjects, and each sample was tested twice by three different lots (total of six results, two replicate tests from one dilution). Aliquots of the venous samples were also tested by the NGSP-certified method, providing approximately 300 comparative results. Data analysis again consisted of least squares linear regression (x = reference results), bias calculation, and Bland Altman limits. The data are provided below.

n	299	Bias at 6% A1C (% difference )	5.77 (- 3.83%)
Slope	1.023	Bias at 7% A1C (% difference)	6.80 (-2.86%)
y-intercept	-0.361	Bias at 9% A1C (% difference)	8.84 (- 1.78%)
"r"	0.90	Avg. % diff.	- 2.82%

The results showed that the accuracy with venous sampling was, on average, 97%. This means that, on average, a true 7 %A1C could read approximately 6.8 %A1C.

An individual result may differ by -1.3 %A1C and +0.9% A1C from the true result. The A1cNow+ system may be used with either fingerstick (capillary) or venous (EDTA anticoagulated) whole blood samples.

## Expected Performance in Waived Laboratories

Clinical studies were performed at three US sites with over 180 untrained people(most with diabetes). These study subjects read the instructions and then performed one A1cNow+ test on themselves. A venous blood sample was collected from each subject, and this sample was tested by an NGSP-certified laboratory method for %A1C. The two results were then compared.

## Untrained User A1cNow+ and an NGSP-certified method

(Tosoh A1c 2.2 Plus)

n	188	Bias at 6% A1C (% difference)	6.02 (+ 0.33%)
Slope	0.99	Bias at 7% A1C (% difference)	7.01 (+0.14%)
y-intercept	0.08	Bias at 9% A1C (% difference)	8.99 (- 0.11%)
"r"	0.93	Avg. % diff.	+ 0.12%

The results showed that untrained users could perform A1cNow+ test-ing on themselves with the same accuracy as trained individuals.

## References

1. Buris, C.A., Ashwood, E.R., Tietz Textbook of Clinical Chemistry, 3rd Edition, W.B. Saunders Co., 1999.

2. Nathan, D.M., et al. The clinical information value of the glycosylated hemoglobin assay. N Engl J Med 1984; 310: 341-346.

3. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complication in insulin-dependent diabetes mellitus. N Engl J Med 1993; 329: 977-986.

4. Diabetes Care 1999; 22 (Suppl. 1); S32-S41.

5. Fogh-Anderson, N., D'Orazio, P.; Proposal for standardizing direct-reading biosensors for blood glucose. Clin Chem 1998; 44(3); 655-659.

6. MLO Supplement, Point-of-Care Testing, 1992.

7. Cagliero, E., Levina, E.V., Nathan, D.M.; immediate feedback of A1C levels improves glycemic control in type 1 and insulin-treated type 2 diabetic patients. Diabetes Care 1999; 22(11): 1785-1789.

8. Goldstein, D.E., Little, R.R., Wiedmeyer, H.M., et al.; Glycated hemoglobin: Methodologies and clinical applications. Clin Chem 1986; 32: B64-B70.

9. American Diabetes Association: Clinical Practice Recommendations 2006. Diabetes Care, 2006; 29 (Suppl. 1)

	INTERNATIONAL SYMBOLS
	MANUFACTURER
	CONTAINS SUFFICIENT FOR <n> TESTS
	IN VITRO DIAGNOSTIC MEDICAL DEVICE
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

