

PRODUCT CATALOGUE



POINT OF CARE TESTS



TABLE OF CONTENTS

PRODUCT GROUP	PAGE
POINT OF CARE READERS	3
DIABETES TESTS	4
CARDIAC MARKERS	4
DRUGS OF ABUSE TESTS	5
ALCOHOL TESTS	6
PREGNANCY TESTS	6
ADULTERATION DETECTION	7
CONSUMABLES	8
HOW TO ORDER	9
ORDER FORM	10

POINT OF CARE READERS

DXPRESS™ POINT OF CARE READER

DXpress™



Product Description	Photometric reader which provides automated result interpretation and documentation for the StatusFirst™ cardiac markers, AccuSign® drugs of abuse and BioSign® hCG tests
Result Type	Quantitative or qualitative
Data Input	Via barcode scanner or keypad
Result Output	Built-in thermal printer or connection to a PC
Result Details	Patient ID; date & time; operator ID; product type; lot number; control line(s); test result(s)
Memory Capacity	Up to 1,000 test results
Power Source	Mains or batteries (6 x AA)
Warranty	1 year (parts & labour)
Regulatory Approval	CE marked for professional use

PRODUCT

PRODUCT CODE

DXpress™ Reader - Purchase (with a 1 year warranty for parts and labour)	LSR-2000-P
DXpress™ Reader - Rental (for a minimum of 1 year)	LSR-2000-R
DXpress™ Reader - Extended Warranty (parts and labour) (for a minimum of 1 year and a maximum of 4 years)	LSR-2000-EW

DIABETES TESTS

A1cNOW+® HbA1c TEST (for healthcare professionals)

A1cNOW+

Product Description	Kit for the monitoring of HbA1c (glycated haemoglobin) levels
Device Format	1 monitor/10 test cartridges/10 sample dilution kits
Sample Medium	Blood (fingerstick or venous)
Sample Volume	10 µL
Result Type	Quantitative (DCCT aligned and NGSP certified)
Result Time	5 minutes
Shelf Life	12 months from the date of manufacture
Storage Temperature	2–8°C - can be stored at room temperature (18–28°C) for up to 3 months
Regulatory Approval	CE marked for professional use



PRODUCT

PRODUCT CODE

PACK SIZE

A1cNow+® HbA1c Test (for healthcare professionals)

03015

10 Test Kit

CARDIAC MARKERS

STATUSFIRST™ TROPONIN I TEST

StatusFirst™

Product Description	Rapid on-site <i>in vitro</i> test for the determination of Troponin I
Device Format	Cassette (single-use)
Sample Medium	Whole blood, serum or plasma
Sample Volume	130 µL
Result Type	Quantitative (range 0.05–30 ng/mL)
Result Time	15 minutes
Shelf Life	12 months from the date of manufacture
Storage Temperature	2–8°C
Regulatory Approval	CE marked for professional use



PRODUCT

PRODUCT CODE

PACK SIZE

StatusFirst™ Troponin I Test (used with DXpress™ Reader)

LSC-131

20 Test Kit

DRUGS OF ABUSE TESTS

ACCUSIGN® DRUGS OF ABUSE TESTS

AccuSign® DOA

Product Description	Rapid on-site <i>in vitro</i> tests for the detection of a range of drugs of abuse
Device Format	Cassette (single-use)
Sample Medium	Urine
Result Type	Qualitative (positive/negative)
Result Interpretation	Using DXpress™ Reader or visually
Result Time	3 or 5 minutes
Shelf Life	15 months from the date of manufacture
Storage Temperature	2–30°C
Regulatory Approval	CE marked for professional use



PRODUCT	PRODUCT CODE	PACK SIZE
AccuSign® Cocaine Test (1 test cassette)	DOA-200-10 DOA-200-35	10 Test Kit 35 Test Kit
AccuSign® Cannabis (THC) Test (1 test cassette)	DOA-201-10 DOA-201-35	10 Test Kit 35 Test Kit
AccuSign® Methamphetamine Test (1 test cassette)	DOA-202-10 DOA-202-35	10 Test Kit 35 Test Kit
AccuSign® Opiates Test (1 test cassette)	DOA-203-10 DOA-203-35	10 Test Kit 35 Test Kit
AccuSign® Amphetamine Test (1 test cassette)	DOA-204-10 DOA-204-35	10 Test Kit 35 Test Kit
AccuSign® Phencyclidine Test (1 test cassette)	DOA-205-10 DOA-205-35	10 Test Kit 35 Test Kit
AccuSign® Barbiturates Test (1 test cassette)	DOA-206-10 DOA-206-35	10 Test Kit 35 Test Kit
AccuSign® Benzodiazepines Test (1 test cassette)	DOA-207-10 DOA-207-35	10 Test Kit 35 Test Kit
AccuSign® Methadone Test (1 test cassette)	DOA-208-10 DOA-208-35	10 Test Kit 35 Test Kit
AccuSign® Tricyclic Antidepressants Test (1 test cassette)	DOA-209-10 DOA-209-35	10 Test Kit 35 Test Kit
AccuSign® Buprenorphine Test (1 test cassette)	DOA-212-10 DOA-212-35	10 Test Kit 35 Test Kit
AccuSign® EDDP (Methadone Metabolite) Test (1 test cassette)	DOA-213-10 DOA-213-35	10 Test Kit 35 Test Kit
AccuSign® 6-Acetylmorphine (Heroin Metabolite) Test (1 test cassette)	DOA-214-10 DOA-214-35	10 Test Kit 35 Test Kit
AccuSign® Cannabis/Opiates/Cocaine/Amphetamine Test (4 test cassette)	DOA-240-10 DOA-240-35	10 Test Kit 35 Test Kit
AccuSign® Cannabis/Opiates/Cocaine/Benzodiazepines/Amphetamine Test (5 test cassette)	DOA-252-10 DOA-252-35	10 Test Kit 35 Test Kit
AccuSign® Amphetamine/Opiates/Cocaine/Benzodiazepines/ Methadone Test (5 test cassette)	DOA-253-10 DOA-253-35	10 Test Kit 35 Test Kit
AccuSign® Amphetamine/Opiates/Cocaine/Benzodiazepines/ Methadone/ Buprenorphine Test (6 test cassette)	DOA-263-10 DOA-263-35	10 Test Kit 35 Test Kit
Test made to individual specifications from any combination of the 1 test cassettes (subject to a minimum order of 3,000 panels per annum)	DOA-XXX-10 DOA-XXX-35	10 Test Kit 35 Test Kit

ALCOHOL TESTS

ACCUSIGN® ALCOHOL TEST

AccuSign®
Alcohol

Product Description	Rapid on-site <i>in vitro</i> test for the semi-quantitative approximation of blood alcohol concentration
Device Format	Dipstick (single-use)
Sample Medium	Saliva
Result Type	Semi-quantitative (blood alcohol concentration range of 0.02%-0.3%)
Result Interpretation	Visual
Result Time	2 minutes
Shelf Life	12 months from the date of manufacture
Storage Temperature	2–27°C



PRODUCT

PRODUCT CODE

PACK SIZE

AccuSign® Alcohol Test (1 test dipstick)	DOA-216-25 DOA-216-50	25 Test Kit 50 Test Kit
--	--------------------------	----------------------------

PREGNANCY TESTS

BIOSIGN® hCG Test (for DXpress™ Reader)

BioSign® hCG

Product Description	Rapid on-site <i>in vitro</i> test for the detection of hCG (human chorionic gonadotropin)
Device Format	Cassette (single-use)
Sample Medium	Urine
Sample Volume	3 drops (110 µL)
Result Type	Qualitative (positive/negative) with borderline result (cutoff concentration 25 mIU/mL)
Result Time	5 minutes
Shelf Life	24 months from the date of manufacture
Storage Temperature	2–8°C
Regulatory Approval	CE marked for professional use



BIOSIGN® hCG Test (for visual interpretation)

BioSign® hCG

Product Description	Rapid on-site <i>in vitro</i> test for the detection of hCG (human chorionic gonadotropin)
Device Format	Cassette (single-use)
Sample Medium	Urine
Sample Volume	3 drops (110 µL)
Result Type	Qualitative (positive/negative) result (cutoff concentration 15 mIU/mL)
Result Time	3 minutes
Shelf Life	24 months from the date of manufacture
Storage Temperature	2–8°C
Regulatory Approval	CE marked for professional use



PRODUCT

PRODUCT CODE

PACK SIZE

BioSign® hCG (Pregnancy) Test (for DXpress™ Reader) (1 test cassette)	BSP-121-10R BSP-121-35R	10 Test Kit 35 Test Kit
BioSign® hCG (Pregnancy) Test (for visual interpretation) (1 test cassette)	BSP-121-10V BSP-121-35V	10 Test Kit 35 Test Kit



ADULTERATION DETECTION

INTECT® 7

Product Description	One step in vitro test for the detection of creatinine, nitrite, glutaraldehyde, pH, specific gravity, bleach and pyridinium chlorochromate
Device Format	Dipstick (single-use)
Sample Medium	Urine
Result Type	Semi-quantitative
Result Interpretation	Visual
Result Time	2 minutes
Shelf Life	18 months from the date of manufacture
Storage Temperature	15–30°C
Regulatory Approval	CE marked for professional use

Intect® 7



TEMPERATURE STRIP

Test strip for the temperature measurement of urine samples (degrees Fahrenheit scale)
Self-adhesive strip can be attached to any sample collection container
Semi-quantitative result read from illuminated squares
Result time within 4 minutes
A temperature reading outside the range 90°F-100°F suggests that a sample may have been altered or substituted
Room temperature storage
Also available with urine container (see page 8)



TN2 INFRARED THERMOMETER

Non-contact infrared thermometer with laser alignment for the measurement of surface temperature in degrees Celsius or Fahrenheit scales
Result is obtained by pointing the device at the target and pressing the read button
Quantitative result (temperature range: -33°C to 250°C)
Result time within 4 minutes
A reading outside the range 32°C/90°F - 38°C/100°F suggests that a sample may have been altered or substituted (normal body temperature is approx 37°C/98°F)



PRODUCT

PRODUCT CODE





PACK SIZE

Intect® 7 Creatinine/Nitrite/Glutaraldehyde/pH/Specific Gravity/Bleach/ Pyridinium Chlorochromate Test (7 test dipstick)	INT-25	25 Test Kit
Temperature Strip (self-adhesive/90-100°F)	TPS-25 TPS-50	25 Test Kit 50 Test Kit
TN2 Infrared Thermometer	IFT-TN2	1 Device

CONSUMABLES

PRODUCT	PRODUCT CODE	PACK SIZE
Latex Powdered Gloves	GLV-101 (Small)	100 Gloves
	GLV-102 (Medium)	100 Gloves
	GLV-103 (Large)	100 Gloves
Latex Powder Free Gloves	GLV-201 (Small)	100 Gloves
	GLV-202 (Medium)	100 Gloves
	GLV-203 (Large)	100 Gloves
Vinyl Powder Free Gloves	GLV-401 (Small)	100 Gloves
	GLV-402 (Medium)	100 Gloves
	GLV-403 (Large)	100 Gloves
Urine Container with Screw Cap (non-sterile/150mL)	CON-101	50 Containers
Urine Container with Screw Cap and Temperature Strip	CON-102	50 Containers
Urine Container with Screw Cap (sterile/150mL)	CON-201	50 Containers
Open Gallipot (250mL)	GAP-250	50 Gallipots
Open Gallipot (500mL)	GAP-500	50 Gallipots
AccuSign® Drugs of Abuse Tests Negative Control (3mL)	CTR-101	1 Bottle
AccuSign® Drugs of Abuse Tests Positive Control (3mL) [custom made to any combination of drugs required]	CTR-102	1 Bottle
BioSign® hCG Test Negative Control (3mL)	CTR-201	1 Bottle
BioSign® hCG Test Positive Control (3mL)	CTR-202	1 Bottle
Printer Paper - for DXpress™ Reader (Width 57mm; Diameter 30mm; Length 10.5m)	TPR-12	12 Rolls
	TPR-24	24 Rolls
Batteries - for DXpress™ Reader (AA, 1.5V Alkaline)	BAT-101	12 Batteries
Battery - for Infrared Thermometer (CR2032, 3V Lithium)	BAT-201	1 Battery



HOW TO ORDER			
 TEL	<p>Call the Chirus Sales Department on: (01923) 212744 [UK] +44 (0)1923 212744 [International]</p>		
 FAX	<p>Complete the ORDER FORM overleaf or use official purchase order form and fax to: (01923) 244751 [UK] +44 (0)1923 244751 [International]</p>		
 EMAIL	<p>Email the information requested on the ORDER FORM overleaf to: sales@chirus.com</p>		
 POST	<p>Complete the ORDER FORM overleaf or use official purchase order form and mail to:</p> <table border="0" style="width: 100%;"> <tr> <td style="width: 50%; vertical-align: top;"> <p>United Kingdom orders [postage NOT required]</p> <p>Chirus Limited FREEPOST NAT2113 Watford WD18 8BR</p> </td> <td style="width: 50%; vertical-align: top;"> <p>International orders [postage required]</p> <p>Chirus Limited Park House 15-23 Greenhill Crescent Watford Business Park Watford Hertfordshire WD18 8PH United Kingdom</p> </td> </tr> </table>	<p>United Kingdom orders [postage NOT required]</p> <p>Chirus Limited FREEPOST NAT2113 Watford WD18 8BR</p>	<p>International orders [postage required]</p> <p>Chirus Limited Park House 15-23 Greenhill Crescent Watford Business Park Watford Hertfordshire WD18 8PH United Kingdom</p>
<p>United Kingdom orders [postage NOT required]</p> <p>Chirus Limited FREEPOST NAT2113 Watford WD18 8BR</p>	<p>International orders [postage required]</p> <p>Chirus Limited Park House 15-23 Greenhill Crescent Watford Business Park Watford Hertfordshire WD18 8PH United Kingdom</p>		

Notes

- If you have access to a scanner, the completed ORDER FORM can be scanned and converted to an image and sent via email. Most image file formats as well as PDF are acceptable.
- For product prices, use the POINT OF CARE TESTS PRICE LIST, a copy of which is available upon request. Alternatively, contact Chirus for a quote.
- The carriage charge is levied on a weight basis. Contact Chirus for a quote.
- Orders from NHS trusts must be accompanied by a purchase order no if the invoice is to be sent to the trust's finance department.
- Orders received by 3:30pm are dispatched the same day for delivery the following working day. Orders received after 3:30pm are dispatched the following working day. Contact chirus for same day deliveries.
- All goods are sold subject to the Chirus Standard Conditions of Business, a copy of which is available upon request.

ORDER FORM 

PRODUCT	PRODUCT CODE	QUANTITY	PRICE/PACK (exc VAT)	SUBTOTAL (exc VAT)
			TOTAL (exc VAT)	
			CARRIAGE (exc VAT)	Contact Chirus

INVOICE ADDRESS	
Organisation _____	
For the Attention of _____	Unit/Department _____
Address _____	
City/Town _____	County _____ Post Code _____
Tel No _____	Purchase Order No _____
Tick box if you wish to have the invoice included with the order <input type="checkbox"/>	

DELIVERY ADDRESS	
Organisation _____	
For the Attention of _____	Unit/Department _____
Address _____	
City/Town _____	County _____ Post Code _____
Tel No _____	

Sc or photocopy

← HOW TO ORDER