

Nucleic Acid Testing

REF PCR0401

128 tests

M1600-PCR Tube

All trademarks are properties of their respective owners.

Key to Graphical Symbols Used

LOT

batch code



use by



manufacturer



contains sufficient for <n> tests

IVD

in vitro diagnostic medical device



temperature limitation

REF

catalogue number



consult instructions for use

EC **REP**

authorized representative in the European
Community



date of manufacture

EC **REP**

OBELIS S.A.,
Bd. Général Wahis, 53
1030 Brussels
Belgium

Autobio Labtec Instruments Co., Ltd.
No.199,15th Ave
National Eco & Tech Zone
Zhengzhou
China
450016



For any technical assistance please contact us in English at:

Email: customerservice@autobio.com.cn

Contact your local dealers for all product related questions in your local language

Intended use

M1600-PCR Tube is intended used to carry the sample which will be amplified based on the PCR principle on the Automated Nucleic Acid Purification and Real Time PCR System during the test procedure. This M1600-PCR Tube is for professional use only.

Pack specification

Name	Quantity
M1600-PCR Tube	128 tests

Assay Analyzers on which the consumables can be used


M1600-PCR Tube is intended use as the consumable of Automated Nucleic Acid Purification and Real Time PCR System, it can be used on the AutoMolec 1600.

Storage

Store at 0-30°C.

Measurement Procedure

Refer to the AutoMolec 1600 instrument operating manual.



Quality Control