

Nucleic Acid Testing

REF PCR0101

600 tests

M3000-Tip Tray

All trademarks are properties of their respective owners.

Key to Graphical Symbols Used

LOT

batch code



use by



manufacturer



contains sufficient for <n> tests

IVD

in vitro diagnostic medical device



temperature limitation

REF

catalogue number



consult instructions for use

EC **REP**

authorized representative in the European
Community



date of manufacture

EC **REP**

OBELIS S.A.,
Bd. Général Wahis, 53
1030 Brussels
Belgium

Autobio Labtec Instruments Co., Ltd.
No.199,15th Ave
National Eco & Tech Zone
Zhengzhou
China
450016



For any technical assistance please contact us in English at:

Email: customerservice@autobio.com.cn

Contact your local dealers for all product related questions in your local language

Intended use

M3000-Tip Tray is intended used to transfer the purified sample and add the amplification reagent into the PCR tube which will be amplified on the Automated Nucleic Acid Purification and Real Time PCR System during the test procedure.
This M3000-Tip Tray is for professional use only.

Pack specification

Name	Quantity
M3000-Tip Tray	600 tests

Assay Analyzers on which the consumables can be used

M3000-Tip Tray is intended use as the consumable of Automated Nucleic Acid Purification and Real Time PCR System, it can be used on the AutoMolec 3000.

Storage

Store at 0-30°C.

Measurement Procedure

Refer to the AutoMolec 3000 instrument operating manual.

Quality Control

