

Microbiological Assay

REF MG0102

30 tests

2#Tip

All trademarks are properties of their respective owners.

Key to Graphical Symbols Used			
	batch code		use by
	manufacturer		contains sufficient for <n> tests
	<i>in vitro</i> diagnostic medical device		temperature limitation
	catalogue number		consult instructions for use or Consult Electronic Instructions for Use
	authorized representative in the European Conformity		Country and Date of Manufacture

	<p>OBELIS S.A Bd. Général Wahis, 53 1030 Brussels Belgium</p>
	<p>Autobio Labtec Instruments Co., Ltd. No.199,15th Ave National Eco & Tech Zone Zhengzhou 450016 China</p>



For any technical assistance please contact us in English at:

Email: customerservice@autobio.com.cn

Contact your local dealers for all product related questions in your local language

Intended purpose

2#Tip is intended to carry and distribute the sample which will be tested on the Automated Microorganism Identification and Antimicrobial Susceptibility Testing Analyzer during the test procedure.
2#Tip is for professional use only.

Pack specification

Name	Quantity
2#Tip	30 tests

Assay Analyzers on which the consumables can be used

2#Tip is intended to be used as the consumable of Automated Microorganism Identification and Antimicrobial Susceptibility Testing Analyzer(AutoMic-i600), and it can be used on the analyzer.

Storage

Store at 10-40°C.

Measurement Procedure

Refer to the operating manual of AutoMic-i600.

Quality Control

